



download software from www.mindresearch.xyz

Documentation (release 4.3)


Contents:

1. Download and Installation of Muse_EEG_View
2. Installing "Winsock Listener" on Port 5000
3. Streaming data with MuseMonitor App (smartphone)
4. Streaming data with MuseDirect App (Windows 10)
5. Playback recorded files
6. Live stream Screen
7. Neurofeedback Screen
8. RAW data (voltages)
9. Logarithmic Graph
10. Lineair Graph
11. Relative Graph
12. Asymmetrie Graph
13. Specific waves


1) Installation on Windows PC or Laptop

Download the App


Click the [download button](#) and save the Zip file in any folder you like.



After download: click with right-mouse button and choose to Un-zip the file




In downloads you see a Muse_EEG_V43 folder and beneath that a <Muse-EEG-View> folder.




Copy or Drag this folder to a data-drive like D:\ or E:\ etc.

You can copy the <Muse EEG View> folder to any place on your computer. But please do **not** use C:\ because this drive may refuse installation of the so called WinSocked component that is needed for connection to MuseMonitor or MuseDirect

Start the App click – with right mouse on EXE file



the first time choose to run as administrator



2. Installing “Winsock Listener” on Port 5000

To listen to Port 5000 (where the streaming data comes in) the App needs a so called WinSock component. This component is missing on most Window PC/Laptops

In the installation packet the MSWINSCK.OCX component is included, but you have to tell windows here-I-am. This registration has to be done only once. If you start the live-streaming without registration you will get a Error message.

How to register the OCX

1. Start the App with administrator permissions ! This is done by clicking with right-mouse on the Muse-EEG-View icon (or muse-eeg-view.EXE) and choose <Run as Administrator>
2. Then click the <register MSWINSCK button> on the startup screen.
If you get ErrNum 75 that means you have no administrator-write-permissions!
If registration is oke you probably get the message : Registration OCX succeeded
3. Now you can start the <Live stream> screen.
If the screen appears without a error the OCX registration is oke.
4. If the above does not work and you did a installation on drive C:\ then please <drag or copy> the <Muse_EEG_View> folder to a data-drivre like D:\ or E:"
5. If you start streaming and no data appear there is a diagnostic screen to monitor every " bit and byte" that is received on port-5000.

3. Streaming data with the MuseMonitor App



Windows version 7 - 10

A stream from MuseMonitor works on Windows7-10, MuseDirect is only available for Windows 10.

Make sure you have the new version of MuseMonitor App on your smartphone.
Because in the old version only ONE value was given for each wave (delta, beta etc)
In the new version a value for each wave on each elektrode is transmitted.
If the MuseMonitor App is running go to settings (the little wheel)

To stream data from MuseMonitor I use these settings:

- * Show values = activated
- * Recording interval - while streaming the interval setting does not matter
- * Recording format : avarage
- * OSC Stream IP - the IP-adress of your PC/Laptop
- * OSC Stream IP – default = 5000
- * OSC stream brainwaves : All (not avarage)
- * Graph type = plot absolute waves ...
- * Brain-waves(absolute)
Activate : Delta - Theta – Alpha – Beta - Gamma
- * Menu Wave-type (absolute)
Average values for ALL sensors TP9 AF7AF8 TP10)

- * Advances settings
Check notch -> in USA = 60


Start streaming

There is a button on MuseMonitor and the word “streaming” appers.

4. Streaming data with MuseDirect (on Windows 10)

The MuseDirect App is for free. It is the substitute for Muse-IO used on headband-I (2014).

Add a new port -5000 (any name you like)



Settings (of the new port)

poort5000

OSC UDP

IP: 127.0.0.1 Port: 5000

Prefix

- Person's Name (e.g. Person1)
- Muse ID (e.g. Muse-1234)
- Custom Static Text:

Example: Muse-1D5C/eeg/0

Output Data:

- EEG
- Accelerometer
- Gyroscope
- Battery

Output Algorithm:

- Absolute Band Powers
- Relative Band Powers
- Band Power Scores
- HSI Precision

5) Playback recorded files


You can playback files that are recorded on the “live screen” in Muse-EEG-View.
You can also playback CSV files that were recorded with MuseMonitor.

a. CSV Files from MuseMonitor

Please put the CSV files from muse-monitor in the Muse_EEG_view/data folder.
They will show up in the file-display-box in the App with the extension CSV.
After you click on it – and graph it for the first time - the App will do a file-conversion.
After that it gets a TXT extension and the layout is changed.
The original file is placed in the <save-csv> folder so you can always retrieve it.

b. Possibilities on the Playback screen

Two windows



There are 2 sub-screens in which you can draw a graph and switch between the two.
That makes it easy to compare 2 graphs or persons for example :
** asymmetry on AF7-AF8 in the upper versus TP9-TP10 in lower screen
** the relative eef waves of Maria versus those of James
** compare high peaks in the delta-wave with bad-disconnection in the raw-data graph

Quick display / Playback

With quick display the graph is drawn directly and stops at the end of the screen. You have to click <continue> to go to the next screen.



With “play back” activated the recording will be played back in the same ‘time dimension’ and you may click <pause> to halt for a moment

Graph Width / Amplitude



With some people the activity is quite strong and gets out of the picture-box. Then you can make the amplitude a bit smaller. During meditation the band-power (waves) are smaller and you might want to amplify a bit to see more distinction.

With a recording of 256Hz (data-points per second) the separate points are drawn very close together and you might give them more stretch/width. Or you want to see the result of 2 screens on 1 screen then set the width to 2 and the data-points are drawn closer together.

From time – until time

From minute/second	Until minute/second
0 0	60 0

You can choose a certain time-interval. The graph begins at starting-time and stops at end-time. Also the (average) calculations are done only over this period. This comes in handy when you do a recording with different activity's and what to see the changes between them. Or when there is a bad connection over a period of time which you want to exclude in the calculations.

Compress files

Compress File ->	Compression is handy if you want to see the average display of a long session. You can choose the compression-time pro 1 seconds or 10 or 30 or 1 minuut
------------------	--

For example you did a recording of half hour with measurement every 0.5 second and want to see the values pro minuut on one screen. Or you did a <constant> recording of 256 measurements pro second and want to see the results pro second.

The average is calculated and written in a new file. The original file will is not changed To view/select the compression click on <show compressed files>.

When the name begint with C1- it is a second-compression, C10 is 10 second compression etc So you can do more time-compressions on one file c.q. one recording session.


6. Live stream Screen en Neurofeedback

(Above: snapshot of neurofeedback screen)

The possibility's of "Live EEG stream" are the same as on the playback screen
You can do neurofeedback on a specific wave-band (delta to gamma) or on the synchronicity (symmetry) between eleft-and-right hemisphere.
Feedback can be done by moving a boat on a lake. Other possibility's are coming soon ...

8) RAW voltages


This graph shows the so called RAW voltages (mVolt). They may vary between 600 and 1000mV



Bad connection


During recording there may be a bad connection on 1 of the 4 electrodes. In the graph you may see much up-and-downs (heavy swing). With a bad connection the App changes the draw-thickness of the graph-line of that electrode and make it thinner so you may recognize the bad connection (in the raw-voltage graph).

9) Logaritmic Graph




The power of a frequency range (for instance alpha 9-13Hz) is called absolute linear band power. First these are calculated in units dB (f.i. the average alpha-power is 60 dB). After that – to make a nicer graph – the Logaritm values are calculated. In the CSV file the band-power is in logaritm values. Some of these will be negative (i.e. when the absolute linear power is less than 1).

10) Linear Graph



To draw this graph the logarithm values are transformed to their 'original' linear values. The linear-scale band power can be calculated from the log-scale band power thusly:
 $\rightarrow \text{linear-scale band power} = 10^{\text{(log-scale band power)}}$.


11) Relative Graph



Researchers often use the relative power (percentages) as a measurement in their research. These are calculated as percentages of linear band powers in each band. The relative-values is the power of 1 band (wave) compared to the total of all bands. i.e. $\alpha\text{-relative-\%} = \left(\frac{\alpha}{\text{delta} + \text{theta} + \alpha + \text{beta} + \text{gamma}} \right) * 100$

12) Asymmetrie Graph

When our state-of-mind changes by the actions we perform not only the power (amplitude) of the waves change but also the contribution of right and left part of the brain (hemisphere). A-symmetry is the (relative) difference in power between left and right hemisphere of the brain. These data are important when researching i.e. meditation, concentration, learning problems etc. Want to know more? Search on Google for <alpha power symmetry index>



	Difference Linear values	Symmetrie Log(..) * 10
TP10 TP9	95%	-0,67
TP10 TP9	100%	-0,12
TP10 TP9	98%	-0,26
TP10 TP9	87%	-1,35
TP10 TP9	68%	-3,58

There are two graphs in the App that show the asymmetrie in power of different waves.

- * a graph for the frontal channels AF8 (right) compared to AF7 (left)
- * a graph for the pariental channels TP10 (right) compared to TP9 (left)

There are several algorithms (ways of calculating) these differences.

The most common is: take the logarithm of right divided by left

for example: index-frontal-alpha-asymmetrie = $\log(\text{alpha-AF8}/\text{alpha-AF7})$

	Difference Linear values	Symmetrie Log(..) * 10	
AF8 AF7	74%	-3,67	<p>The percentage of correspondance is the % of difference in the total linear power right/left. i.e. if beta on TP9 is 50dB and on TP10 is 45dB their correspondance is 90%</p> <p>The average symmetrie over the session.</p> <ul style="list-style-type: none"> • If this logarithm has a positive (+) value it means more activity at the right site. • If the logarithm has a negative (-) value it means more activity at the left site.
AF8 AF7	72%	-3,70	
AF8 AF7	91%	0,80	
AF8 AF7	94%	1,19	
AF8 AF7	88%	1,74	
AF8 AF7			

13) Specific waves

Specific waves

Select Data


- All Power Bands
- All Power Bands
- TP9 (5 powerbands)
- AF7 (5 powerbands)
- AF8 (5 powerbands)
- TP10 (5 powerbands)
- Delta (4 electrodes)
- Theta (4 electrodes)
- Alpha (4 electrodes)

If you want to dive deeper into a specific problem you can:

- * display a single waves on 4 different electrodes
- * display 5 bands/waves on 1 single electrode/channel

Suppose you want to know where the high-alpha comes from ?
 Select alpha in the < select data > box and you see alpha on 4 the electrodes. You may find that delta on AF7 is very high (due to a bad connection of the haedband)

Why is Delta so high ?



The Muse headset has relative high Delta values. As James - the developer of MuseMonitor - explains there are 2 reasons :

a) During eye movement there are electrical (non brain) signals that make big slow (relatively) arcs in the raw data. These electrical currents interfere with the brain data, especially in the lowest frequency band and causes Delta to rise. So – while recording a session – try to avoid eye movements.

b) It is due to the kind of FFT (fast fourier) calculation the Muse headset uses. Because the EEG wave bands represent a power spectrum, their values will vary exponentially, meaning the lower-frequency bands (such as delta and theta) will be exponentially larger values than the higher-frequency bands (alpha and beta).

<< End >>